


## Terragene Sterilisation Control - Biological Indicators

Self-Contained biological indicators for monitoring steam sterilisation processes



## Features

Biological indicators are the only internationally accepted devices that provide a direct measure of a sterilisation cycle lethality.

### Self-Contained Biological Indicators

- With ever increasing demands for faster sterile instruments turnaround in healthcare facilities, sterilisation results need to be available ASAP to reliably verify the sterility of the load.
- Terragene Bionova Ultra Rapid and Super Rapid Biological Indicators for Steam, Ethylene Oxide, Hydrogen Peroxide and Formaldehyde sterilisation processes provides you with reliable results in record time, to ultimately speed up your workflow.
- Terragene Bionova Ultra Rapid and Super Rapid innovative technology is one of the fastest biological indicators on the market.
- Terragene Bionova biological indicators provide you with final readouts within 20 minutes to 4 hours, depending on type.

### Conventional Biological Indicators

- Using Terragene Bionova Conventional SCBI's in conjunction with Terragene's Bionova Dual Incubator, you can produce viable results with 24 h for Steam & VH<sub>2</sub>O<sub>2</sub> or 48 h for Ethylene Oxide, Formaldehyde and Dry Heat.

Product Code	Product Description	Pack Size	Quote
TKT-BT224	Terragene Bionova Ultra Rapid, Steam, 20 mins	10 Pk or 50 Pk	<a href="#">Download</a>
TKT-BT222	Terragene Bionova Super Rapid, Steam, 1 hr	50 Pk	<a href="#">Download</a>
TKT-BT20	Terragene Bionova, Steam, 24 hr	100 Pk	<a href="#">Download</a>
TKT-BT96	Terragene Bionova Super Rapid, VH <sub>2</sub> O <sub>2</sub> , 30 mins	10 Pk or 50 Pk	<a href="#">Download</a>
TKT-BT95	Terragene Bionova, VH <sub>2</sub> O <sub>2</sub> , 2 hr	50 Pk	<a href="#">Download</a>
TKT-BT91	Terragene Bionova, VH <sub>2</sub> O <sub>2</sub> , 24 hr	50 Pk	<a href="#">Download</a>
TKT-BT30	Terragene Bionova, Dry Heat, 48 hr	50 Pk	<a href="#">Download</a>
TKT-BT110	Terragene Bionova Rapid, Ethylene Oxide, 4 hr	50 Pk	<a href="#">Download</a>
TKT-BT10	Terragene Bionova, Ethylene Oxide, 48 hr	30 Pk or 100 Pk	<a href="#">Download</a>
TKT-BT102	Terragene Bionova Rapid, Formaldehyde, 2 hr	50 Pk	<a href="#">Download</a>
TKT-BT100	Terragene Bionova, Formaldehyde, 48 h	100 Pk	<a href="#">Download</a>

## Relevant Product

---

- [Washing Control - Washer Disinfectors](#)
- [Washing Control - Ultrasonic Cleaning Performance](#)
- [Washing Control - Ultrasonic Cavitation Performance](#)
- [Disinfection Control - Thermal](#)
- [Hygiene Control - WD Protein Detection](#)
- [Sterilisation Control - Chemical Indicators](#)
- [Hygiene Control - Detection Systems](#)
- [Smeg WD2145](#)
- [Smeg WD1160](#)
- [Smeg WD4060](#)
- [Smeg WD4190](#)
- [UW565 Ultrasonic Cleaner](#)
- [Bandelin Sonorex](#)
- [Bandelin Sonomic](#)
- [Bandelin Triston](#)
- [Bandelin Sonoboard](#)

## Applications

---

- Dental
- CSD
- Healthcare Facilities

## Resources

---

- Video
- Video
- Video
- Video
- Video
- Video
- Brochure